

# GAMT antibody



Catalog Number: 107288

## Product name

GAMT antibody

## Specificity

Human, Rat, Mouse, Pig; other species not tested.

## Antibody description

GAMT Mouse Monoclonal antibody. Positive WB detected in fetal human brain tissue, mouse liver tissue, pig liver tissue, rat liver tissue. Observed molecular weight by Western-blot: 26 kDa

## Preparation

This antibody was obtained by immunization of GAMT recombinant protein (Accession Number: NM\_000156). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107288(GAMT Antibody) at dilution of 1:2000