

Human PCSK9/NARC1/PSEC0052 (His/HA/AVI tag) recombinant protein



Catalog Number: 519288

General Information

Gene Name Synonym

Proprotein Convertase Subtilisin/Kexin Type 9, Neural Apoptosis-Regulated Convertase 1, NARC-1, Proprotein Convertase 9, PC9, Subtilisin/Kexin-Like Protease PC9, PCSK9, NARC1

Protein Construction

A DNA sequence encoding the Human PCSK9/NARC1/PSEC0052 Gln31-Gln692 is expressed with a 6His/HA/AVI tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 50mM HEPES, 150mM NaCl, pH 7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PCSK9/NARC1/PSEC0052 (His/HA/AVI tag) recombinant protein