

Human SMR3B (Fc Tag) recombinant protein



Catalog Number: 500688

General Information

Gene Name Synonym

Proline-rich peptide P-B; Proline-rich protein 3;
Peptide P-A; Peptide D1A

Protein Construction

A DNA sequence encoding the human SMR3B (NP_006676.1) (Met1-Pro79) was expressed with the Fc region of mouse IgG1 at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 23

Molecular Mass

The recombinant human SMR3B/mFc comprises 291 amino acids and has a predicted molecular mass of 32.2 kDa. The apparent molecular mass of the monomer is approximately 37 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human SMR3B Protein (Fc Tag) SDS-PAGE