

FLJ14213 antibody



Catalog Number: 110688

Product name

FLJ14213 antibody

Specificity

Human; other species not tested.

Antibody description

FLJ14213 Rabbit Polyclonal antibody. Positive WB detected in DU 145 cells. Observed molecular weight by Western-blot: 27 kDa

Preparation

This antibody was obtained by immunization of FLJ14213 recombinant protein (Accession Number: BC008922). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



DU 145 cells were subjected to SDS PAGE followed by western blot with Catalog No:110688(FLJ14213 Antibody) at dilution of 1:300