

TSPAN7 Polyclonal Antibody



Catalog Number: 161688

Product name

TSPAN7 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to TSPAN7

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 110-249 of human TSPAN7 (NP_004606.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - TSPAN7 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using TSPAN7 antibody . Blue: DAPI for nuclear staining.