

Anti-Glypican 3/GPC3/OCI-5 antibody



Catalog Number: 104688

Product name

Monoclonal

Anti-Glypican 3/GPC3/OCI-5 antibody

Ig Type

Specificity

Rabbit IgG

Human Glypican 3 / GPC3 / OCI-5

Applications

Antibody description

FCM

Rabbit monoclonal to Glypican 3/GPC3/OCI-5

Dilutions

FCM: 0.5-2 µg/Test

Preparation

Validations

This antibody was obtained from a rabbit immunized with a synthetic peptide corresponding to the center region of the Human Glypican 3 / GPC3 / OCI-5, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS


Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality



Glypican 3 / GPC3 / OCI-5 Antibody, Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Human GPC3 expression on HepG2 cells. Cells were stained with purified anti-Human GPC3, then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.