

# C1QTNF3 antibody

Catalog Number: 108688

## Product name

C1QTNF3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

C1QTNF3 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human small intestine tissue.

## Preparation

This antibody was obtained by immunization of C1QTNF3 recombinant protein (Accession Number: BC016021). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:108688(C1QTNF3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:108688(C1QTNF3 Antibody) at dilution of 1:50 (under 40x lens)