

Mouse Transferrin Receptor/TFRC/CD71 (His Tag) recombinant protein



Catalog Number: 501388

General Information

Protein Construction

A DNA sequence encoding the extracellular domain of mouse TFRC (Q62351) (Cys 89-Phe 763) was fused with a polyhistidine tag at the N-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 92 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant mouse TFRC comprises 694 amino acids and predicts a molecular mass of 78.3 kDa. The apparent molecular mass of mouse TFRC

is approximately 80 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse Transferrin Receptor / TFRC / CD71 Protein (His Tag) SDS-PAGE