

PSMC1 antibody



Catalog Number: 114388

Product name

PSMC1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PSMC1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, A549 cells, HeLa cells, Jurkat cells, mouse lung tissue. Observed molecular weight by Western-blot: 52-57 kDa

Preparation

This antibody was obtained by immunization of PSMC1 recombinant protein (Accession Number: NM_002802). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114388(PSMC1 Antibody) at dilution of 1:1000