

# MTMR6 antibody



Catalog Number: 112788

## Product name

MTMR6 antibody

IHC: 1:20-1:200


## Specificity

Human, Mouse, Rat; other species not tested.

## Validations

## Antibody description

MTMR6 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, mouse heart tissue. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 64-67 kDa, 72-76 kDa




## Preparation

This antibody was obtained by immunization of MTMR6 recombinant protein (Accession Number: NM\_004685). Purification method: Antigen affinity purified.

human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112788(MTMR6 antibody) at dilution of 1:500

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemical of paraffin-embedded human liver using Catalog No:112788(MTMR6 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:112788(MTMR6 antibody) at dilution of 1:100 (under 40x lens)