

# PGLYRP4 antibody



Catalog Number: 113788

## Product name

PGLYRP4 antibody

## Specificity

Human; other species not tested.

## Antibody description

PGLYRP4 Rabbit Polyclonal antibody. Positive IHC detected in human skin tissue, human oesophagus tissue.

## Preparation

This antibody was obtained by immunization of PGLYRP4 recombinant protein (Accession Number: NM\_020393). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:113788(PGLYRP4 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:113788(PGLYRP4 Antibody) at dilution of 1:50 (under 10x lens)