

# Anti-SQSTM1/p62 Rabbit antibody



Catalog Number: 167788

## Product name

Anti-SQSTM1/p62 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

SQSTM1/p62 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human SQSTM1/p62

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


WB: 1/1000

ICC/IF: 1/20

FC: 1/50

IP: 1/20

## Validations



Western blot detection of SQSTM1/p62 in A549,HL-60,U2OS cell lysates using SQSTM1/p62 Rabbit mAb(1:1000 diluted).Predicted band size:48kDa.Observed band size:62kDa.