

ADH4 antibody

Catalog Number: 107788

Product name

ADH4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ADH4 Rabbit Polyclonal antibody. Positive IHC detected in human kidney. Positive WB detected in mouse liver tissue, HeLa cells, mouse brain tissue. Observed molecular weight by Western-blot: 43kd

Preparation

This antibody was obtained by immunization of ADH4 recombinant protein (Accession Number: NM_000670). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107788(ADH4 antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human kidney using Catalog No:107788(ADH4 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:107788(ADH4 antibody) at dilution of 1:50 (under 40x lens)