

NEDD8 antibody

Catalog Number: 113098

Product name

NEDD8 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NEDD8 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in mouse brain tissue, human heart tissue. Observed molecular weight by Western-blot: 6-10 kDa

Preparation

This antibody was obtained by immunization of NEDD8 recombinant protein (Accession Number: NM_006156). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113098(NEDD8 antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:113098(NEDD8 Antibody) at dilution of 1:50 (under 40x lens)