

Anti-AMPK alpha 1/2 Rabbit antibody



Catalog Number: 167098

Product name

Anti-AMPK alpha 1/2 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, IP

Antibody description

AMPK alpha 1/2 Rabbit polyclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Recombinant protein of human AMPK alpha 2

IP: 1/20

Validations

Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal



Western blot detection of AMPK alpha 1/2 in Jurkat cell lysates using AMPK alpha 1/2 Rabbit pAb(1:500 diluted). Predicted band size:64,62kDa.Observed band size:62kDa.