

ACPP antibody

Catalog Number: 107098

Product name

ACPP antibody

Specificity

Human; other species not tested.

Antibody description

ACPP Mouse Monoclonal antibody. Positive WB detected in Human testis. Positive IHC detected in human prostate hyperplasia tissue. Observed molecular weight by Western-blot: 45-50kd

Preparation

This antibody was obtained by immunization of ACPP recombinant protein (Accession Number: NM_001134194). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107558(ACPP antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human prostate hyperplasia using Catalog No:107558(ACPP antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human prostate hyperplasia using Catalog No:107558(ACPP antibody) at dilution of 1:50 (under 40x lens)