

Anti-GST-tag Mouse mAb (HRP Conjugate)



Catalog Number: 169098

Product name

Mouse IgG2b

Anti-GST-tag Mouse mAb (HRP Conjugate)

Applications

Antibody description

WB, ELISA

GST-tag Mouse mAb (HRP Conjugate)

Dilutions

Preparation

WB:1/500—1/5000

Antigen: GST fusion protein

Validations

Formulation

In PBS(pH 7.4) containing with 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of purified recombinant GST fusion proteins with different molecular weights using GST-tag Mouse mAb (HRP Conjugate) (1:1000 diluted).