

# Anti-GST-tag Mouse mAb (HRP Conjugate)



Catalog Number: 169098

## Product name

Mouse IgG2b

Anti-GST-tag Mouse mAb (HRP Conjugate)

## Applications

## Antibody description

WB, ELISA

GST-tag Mouse mAb (HRP Conjugate)

## Dilutions

## Preparation

WB:1/500—1/5000

Antigen: GST fusion protein

## Validations

## Formulation

In PBS(pH 7.4) containing with 50% glycerol.


## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of purified recombinant GST fusion proteins with different molecular weights using GST-tag Mouse mAb (HRP Conjugate) (1:1000 diluted).