

Matrilin 2 antibody



Catalog Number: 112498

Product name

Matrilin 2 antibody

Immunogen

[Human MATN2 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

Matrilin 2 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, HEK-293 cells, mouse liver tissue. Observed molecular weight by Western-blot: 120-150 kDa

Preparation

This antibody was obtained by immunization of Matrilin 2 recombinant protein (Accession Number: NM_030583). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:112498(MATN2 antibody) at dilution of 1:800