

Anti-VEGF Rabbit antibody

Catalog Number: 167498

Product name

Anti-VEGF Rabbit antibody

Specificity

Human

Antibody description

VEGF Rabbit polyclonal antibody

Preparation

Antigen: The antiserum was produced against synthesized peptide derived from the Internal region of human VEGFA. AA range:101-150

Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

-20°C/1 year

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ELISA


Dilutions

WB: 1/500 - 1/2000


1:100-300

ELISA: 1/20000

Validations




Immunohistochemical analysis of paraffin-embedded human-heart, antibody was diluted at



1:100

Immunohistochemical analysis of paraffin-embedded human-heart, antibody was diluted at



1:100

Western Blot analysis of 22RV-1, SW480, BT474, K562, mouse liver cells using VEGF-A Polyclonal Antibody. Antibody was diluted at 1:500. Secondary antibody was diluted at 1:20000