

COX7A2 antibody



Catalog Number: 109498

Product name

COX7A2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

COX7A2 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, rat brain tissue. Positive IF detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 9 kDa

Preparation

This antibody was obtained by immunization of COX7A2 recombinant protein (Accession Number: NM_001865). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109498(COX7A2 Antibody) at dilution of 1:300



Immunofluorescent analysis of (-20oc Ethanol) fixed NIH/3T3 cells using Catalog No:109498(COX7A2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)