

# Mouse FTLS / RSPO2 (Fc Tag) recombinant protein



Catalog Number: 500898

## General Information

### Gene Name Synonym

Cysteine-rich and single thrombospondin domain-containing protein 2; Roof plate-specific spondin-2

### Protein Construction

A DNA sequence encoding the mouse Rspo2 (NP\_766403.1) (Met1-Gly205) was expressed with the Fc region of mouse IgG1 at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

> 85 % as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Asn 24

### Molecular Mass

The recombinant mouse Rspo2 consists 416 amino acids and predicts a molecular mass of 47.4 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at -20°C to -80°C.

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse R-Spondin 2 / FTLS / RSPO2 Protein (Fc Tag) SDS-PAGE