

Human SCG3 / Secretogranin 3 (His Tag) recombinant protein



Catalog Number: 502898

General Information

Gene Name Synonym

Secretogranin III

Protein Construction

A DNA sequence encoding the human SCG3 (AAH14539.1) (Met1-Leu468) was expressed with a C-terminal polyhistidine tag.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Phe 20

Molecular Mass

The recombinant human SCG3 comprises 460

amino acids and has a predicted molecular mass of 52.4 kDa. The apparent molecular mass of the protein is approximately 54-63 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human SCG3 / Secretogranin 3 Protein (His Tag)
SDS-PAGE