

# AEBP2 antibody

Catalog Number: 107898

## Product name

AEBP2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

AEBP2 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 52 kDa

## Preparation

This antibody was obtained by immunization of AEBP2 recombinant protein (Accession Number: BC015624). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107898(AEBP2 antibody) at dilution of 1:100



IP Result of anti-AEBP2 (IP:Catalog No:107898, 4ug; Detection:Catalog No:107898 1:300) with mouse testis tissue lysate 4800ug.