

Desmocollin 2 antibody



Catalog Number: 109898

Product name

Desmocollin 2 antibody

Specificity

Human; other species not tested.

Antibody description

Desmocollin 2 Rabbit Polyclonal antibody. Positive IHC detected in human oesophagus tissue.

Preparation

This antibody was obtained by immunization of Desm^oCollin 2 recombinant protein (Accession Number: NM_024422). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:109898(DSC2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:109898(DSC2 Antibody) at dilution of 1:50 (under 40x lens)