

Human DAPK3/ZIPK (GST Tag) recombinant protein



Catalog Number: 502198

General Information

Protein Construction

A DNA sequence encoding the full length of human DAPK3 (NP_001339.1) (Met 1-Arg 454) was fused with the GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

The specific activity was determined to be 5 nmol/min/mg using MBP as substrate.

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human DAPK3/GST chimera

consists of 678 amino acids and predicts a molecular mass of 79 kDa. It migrates as an approximately 70 kDa band in SDS-PAGE under reducing conditions.

Formulation

Supplied as sterile 20mM Tris, 500mM NaCl, 10mM GSH, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human DAPK3 / ZIPK Protein (GST Tag) SDS-PAGE