

Anti-CSN1S1 antibody



Catalog Number: 175198

Product name

Anti-CSN1S1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CSN1S1

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human alpha casein 31-130/185

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB:1:500-2000

Validations




Fig1: Sample:; Human breast milk Lysate at 2 uL; Human breast milk Lysate at 6 uL; Human breast milk Lysate at 8 uL; Primary: Anti-Alpha casein at 1/2000 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 28 kD; Observed band size: 21 kD