

FAM71F2 antibody



Catalog Number: 110598

Product name

FAM71F2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FAM71F2 Rabbit Polyclonal antibody. Positive IF detected in PC-3 cells. Positive WB detected in PC-3 cells, mouse liver tissue. Observed molecular weight by Western-blot: 40-45 kDa

Preparation

This antibody was obtained by immunization of FAM71F2 recombinant protein (Accession Number: XM_011516132). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110598(FAM71F2 Antibody) at dilution of 1:2000



Immunofluorescent analysis of (-20oc Ethanol) fixed PC-3 cells using Catalog No:110598(FAM71F2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)