

Rhesus CD16 / FCGR3 (His & AVI Tag) recombinant protein



Catalog Number: 504598

General Information

Protein Construction

A DNA sequence encoding the rhesus CD16 (NP_001258584.1) extracellular domain (Met 1-Gln 208) was fused with a c-terminal polyhistidine tagged AVI tag at the C-terminus.

Organism

Rhesus

Expression Host

Human Cells

QC Testing

Activity

Measured by its binding ability in a functional ELISA.

Immobilized Rhesus CD16-AVI-His at 10 µg/ml (100 µl/well) can bind biotinylated human IgG1. The EC₅₀ of biotinylated human IgG1 is 0.2-0.5 µg/ml.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gly 17

Molecular Mass

The recombinant rhesus CD16 consists of 223 amino acids and has a calculated molecular mass of 25.3 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Cynomolgus CD16 / FCGR3 Protein (His & AVI Tag) SDS-PAGE