

HTT Polyclonal Antibody



Catalog Number: 165598

Product name

HTT Polyclonal Antibody

Specificity

Rat

Antibody description

Polyclonal antibody to HTT

Preparation

Antigen: A synthetic peptide of human HTT

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

IHC

Dilutions

IHC 1:50 - 1:200

Validations

Immunohistochemistry - HTT Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using HTT antibody at dilution of 1:100 (40x lens).