

# Anti-UBASH3A antibody



Catalog Number: 106598

## Product name

Anti-UBASH3A antibody

## Immunogen

[Human UBASH3A \(aa 354-623, His Tag\) recombinant protein](#)

## Specificity

Human UBASH3A

## Antibody description

Rabbit polyclonal to UBASH3A

## Preparation

Produced in rabbits immunized with purified, recombinant Human UBASH3A (rh UBASH3A; AAH69511.1; Ala354-Asn623). UBASH3A specific IgG was purified by Human UBASH3A affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC-P


## Dilutions

ELISA:0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human UBASH3A. The detection limit for Human UBASH3A is < 0.039 ng/well.

IHC-P: 0.1-2 µg/ml

## Validations



UBASH3A Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry