

Human PFKFB1/F6PK/PFRX (His tag) recombinant protein



Catalog Number: 518598

General Information

Gene Name Synonym

6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 1, 6PF-2-K/Fru-2,6-P2ase liver isozyme, Fructose-2,6-bisphosphatase, PFKFB1, F6PK, PFRX

Protein Construction

A DNA sequence encoding the Human PFKFB1/F6PK/PFRX Ser2-Tyr471 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PFKFB1/F6PK/PFRX (His tag)
recombinant protein