

GNB4 antibody

Catalog Number: 110998

Product name

GNB4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GNB4 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse brain tissue, SH-SY5Y cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 38 kDa

Preparation

This antibody was obtained by immunization of GNB4 recombinant protein (Accession Number: XM_047448653). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110998(GNB4 antibody) at dilution of 1:500



IP Result of anti-GNB4 (IP:Catalog No:110998, 4ug; Detection:Catalog No:110998 1:800) with mouse brain tissue lysate 4000ug.