

# Anti-GAB1 Rabbit antibody

Catalog Number: 166998

## Product name

Anti-GAB1 Rabbit antibody

## Specificity

Human, Rat

## Antibody description

GAB1 Rabbit polyclonal antibody

## Preparation

Antigen: Peptide sequence around aa.625~629 (V-E-Y-L-D) derived from Human Gab1.

## Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB


## Dilutions

WB:1/1000

## Validations



Western blot analysis of extracts from A431 and MCF7 cells using Gab1 (Ab-627) rabbit pAb (1:1000 diluted). Predicted band size:110KDa. Observed band size:110KDa.



Western blot analysis of extracts from A431 and MCF cells using Gab1 (Ab-627) rabbit pAb (1:1000 diluted). Predicted band size:110KDa. Observed band size:110KDa.