

BAG1L-specific antibody



Catalog Number: 106998

Product name

BAG1L-specific antibody

Specificity

Human; other species not tested.

Antibody description

BAG1L-specific Rabbit Polyclonal antibody.
Positive IF detected in HepG2 cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004323).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:10-1:100

Validations



Immunofluorescent analysis of HepG2 cells, using BAG1 antibody Catalog No:108407 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).