

Anti-Matrin 3 Rabbit antibody



Catalog Number: 168998

Product name

Anti-Matrin 3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Matrin 3 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Matrin 3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


1:2000-1:10000

1:200

1:50-1:200

1:50

Validations



Western blot detection of Matrin 3 in Jurkat, HeLa cell lysates using Matrin 3 Rabbit mAb(1:1000 diluted). Predicted band size: 95kDa. Observed band size: 125kDa.