

HERP antibody

Catalog Number: 111298

Product name

HERP antibody

Specificity

Human; other species not tested.

Antibody description

HERP Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells, SW 1990 cells. Positive IHC detected in human lung cancer tissue, human heart tissue, human kidney tissue, human liver tissue, human ovary tissue, human skin tissue, human testis tissue. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of HERP recombinant protein (Accession Number: NM_014685). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:111298(HERPUD1 antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:111298(HERPUD1 antibody) at dilution of 1:100 (under 10x lens)