

NOL9 antibody



Catalog Number: 113298

Product name

NOL9 antibody

Specificity

Human; other species not tested.

Antibody description

NOL9 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells, PC-3 cells. Observed molecular weight by Western-blot: 79 kDa

Preparation

This antibody was obtained by immunization of NOL9 recombinant protein (Accession Number: BC009257). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113298(NOL9 antibody) at dilution of 1:1000