

Anti-S100A7L2 antibody



Catalog Number: 176298

Product name

Anti-S100A7L2 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to S100A7L2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human s100a7l2. aa range:38-87

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of HuvEc cells using S-100A7L2 Polyclonal Antibody