

FOXH1-Specific antibody



Catalog Number: 110698

Product name

FOXH1-Specific antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

FOXH1-Specific Rabbit Polyclonal antibody.
Positive WB detected in mouse liver tissue.
Observed molecular weight by Western-blot:
40-44kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003923).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110698(FOXH1-Specific antibody) at dilution of 1:300