

IGF2BP1 Polyclonal Antibody



Catalog Number: 161698

Product name

IGF2BP1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to IGF2BP1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 480 to the C-terminus of human IGF2BP1 (NP_006537.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF, IP


Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200


IP 1:20 - 1:50

Validations



Western blot - IGF2BP1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using IGF2BP1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunoprecipitation - IGF2BP1 Polyclonal Antibody

Immunoprecipitation analysis of 200ug extracts of K562 cells using 1ug IGF2BP1 antibody. Western blot was performed from the immunoprecipitate using IGF2BP1 antibody at a dilution of 1:1000.