

H3N2 HA (His Tag) recombinant protein



Catalog Number: 503698

General Information

Gene Name Synonym

Hemagglutinin HA1 chain; Hemagglutinin HA2 chain

Protein Construction

A DNA sequence encoding the N-terminal segment (Met 1-Arg 345) of the influenza A hemagglutinin (A/Aichi/2/1968 (H3N2)) (AAA43178.1), termed as HA1, was expressed, fused with a C-terminal polyhistidine tag.

Organism

H3N2

Expression Host

Human Cells

QC Testing

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 17

Molecular Mass

The secreted recombinant HA1 subunit of

influenza A H3N2 (A/Aichi/2/1968 (H3N2)) HA comprises 340 amino acids and has a predicted molecular mass of 37.6 kDa. As a result of glycosylation, it migrates as an approximately 50-55 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H3N2 (A/Aichi/2/1968) Hemagglutinin Protein (HA1 Subunit) (His Tag) SDS-PAGE