

Anti-Bcl-XL Rabbit antibody



Catalog Number: 166698

Product name

Anti-Bcl-XL Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Bcl-XL Rabbit polyclonal antibody

Preparation

Antigen: Peptide sequence around aa.60~64 (A-D-S-P-A) derived from human BCL-XL .

Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB


Dilutions

WB:1/300Other:1/1000

Validations



Western blot detection of BCL2L1 in 3T3,PC12 and HeLa cell lysates using BCL2L1 rabbit pAb (1:1000 diluted).Predicted band size:30KDa.Observed band size:30KDa.



Western blot detection of BCL2L1 in 3T3,PC12 and HeLa cell lysates using BCL2L1 rabbit pAb (dilution 1:300).Predicted band size:30KDa.Observed band size:30KDa.