

UTP6 antibody



Catalog Number: 116698

Product name

UTP6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UTP6 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of UTP6 recombinant protein (Accession Number: NM_018428). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116698(UTP6 Antibody) at dilution of 1:300