

# ACSL5 antibody

Catalog Number: 107698

## Product name

ACSL5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ACSL5 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HeLa cells, mouse liver tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 66-70kd

## Preparation

This antibody was obtained by immunization of ACSL5 recombinant protein (Accession Number: NM\_016234). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107698(ACSL5 antibody) at dilution of 1:1200



IP Result of anti-ACSL5 (IP:Catalog No:107698, 4ug; Detection:Catalog No:107698 1:600) with mouse brain tissue lysate 4000ug.