

NCOA2 Polyclonal Antibody



Catalog Number: 160398

Product name

NCOA2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to NCOA2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 550-820 of human NCOA2 (NP_006531.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - NCOA2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NCOA2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.




Western blot - NCOA2 Polyclonal Antibody

Western blot analysis of extracts from normal (control) and NCOA2 knockout 293T cells, using NCOA2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3min. CRISPR KO cell lines generated by

NCOA2 Polyclonal Antibody



Catalog Number: 160398



Immunohistochemistry - NCOA2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using NCOA2 antibody at dilution of 1:100 (40x lens).