

HECTD3 antibody



Catalog Number: 111398

Product name

HECTD3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

HECTD3 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells. Observed molecular weight by Western-blot: 97kd

Preparation

This antibody was obtained by immunization of HECTD3 recombinant protein (Accession Number: BC019105). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:111398(HECTD3 antibody) at dilution of 1:400