

# TRMT5 antibody



Catalog Number: 116398

---

## Product name

TRMT5 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

human, mouse, rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

TRMT5 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse testis tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 50-55kd

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of TRMT5 recombinant protein (Accession Number: NM\_020810). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116398(TRMT5 antibody) at dilution of 1:500

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality