

IL33 antibody

Catalog Number: 107398

Product name

IL33 antibody

Specificity

Human; other species not tested.

Antibody description

IL33 Mouse Monoclonal antibody. Positive WB detected in human lung tissue, HepG2 cells. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of IL33 recombinant protein (Accession Number: NM_001314045). Purification method: Thiophilic affinity chromatograph.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgM

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:107398(IL33 Antibody) at dilution of 1:1000