

Anti-GAD65 (6D5) Mouse antibody



Catalog Number: 168398

Product name

Anti-GAD65 (6D5) Mouse antibody

Specificity

Human, Rat

Antibody description

GAD65 (6D5) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human GAD2 protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of GAD2 in Rat Hippocampus and Rat Brain lysates using GAD2 mouse mAb(dilution 1:1000).Predicted band size:60kDa.Observed band size:60kDa.